

Dr. Fritz®



GmbH

Endoscopes and Video Systems

Vet-Info 12-3 E

Dogs, Cats

Flexible gastroscopes for examination of the GI tract




Gastroscopy for small animals

Video-Bronchoscopes and Gastrosopes



Completely mounted

**FVE-130-10**

Video gastroscope, 1/6" CCD chip, OD:9,8 mm, WL:130 cm, full size screen technology, standard buttons, freeze function, angle of view 115°, working channel 2,8 mm, standard connectors for cleaning and disinfection, light connector, compatible to Olympus, including density tester, cleaning brush, processing cable, case.

4 way
deflection

Video Processor Light source

– perfect illumination for gastroscopy, laparoscopy and arthroscopy –

Vet Video Combo Xe 100

– the combination of processing unit and Xenon light source universal for all endoscopes –



48-PXE-100



- easy to use, for all applications, economical

48-PXE-100 Vet Video Combo Xe 100, digital CCU & 100 w Xenon light source, for FVE video gastro/colonoscopes and rigid video cameras, PAL, 110-240VAC, 50/60Hz

- light connector comp. to Olympus
- integrated air pump
- automatic exposure control, 2mill.shutter
- 2,5 Lux min, light sensitivity resolution 752(H)x582(V), 440 000 pixel video outlet: VHS, S-VHS,

Accessories:

48-P28-K Camera head

48-XE-100E Xenon spare lamp type XBO R

Suction pump

Suction pump, portable



H160-AB Suction pump, portable complete with 1 Li. bottle (polysulfon) suction fully adjustable vaccum: 0-0,8bar, max. 22 litre (min) 230VAC 50/60Hz

- ideal for the small animal practice
- flexible and easy to use

WV-53022 Video cart, size:62x63cm, h:110cm, completely mounted, basic version, color:grey, partly blue,

- drawer inbuilt
- specially rust protected (without holder for gastroscope and infusion)



Accessories:



48-PK.603 Connecting cord, 18 pin, Dr-Fritz CCU to Dr.Fritz video gastro/colonoscopes, 40cm



X-10-27 Irrigation bottle for gastroscopes comp. to XION and Olympus, 250ml, autoclavable



EF-18 Density tester comp.to flexible video gastro / colonoscopes, Typ FVE



RH99-56 Cap for instrument channel, bag with 5 pcs



F28-1709 Cleaning brush, flexible, 160cm for working channels 2,6- 3,5mm
F28-2208 same but with 220cm










PC-95-PSK USB kit, software driver with cable connection, Chinch/Y/C-USB, comp. to 48-PXE-100







Accessories:

Accessories for high frequency surgery

(Fragen Sie uns nach gebrauchten, leistungsstarken HF-Geräten)

	HF22-1601	Suction-koagulation elektrode , flexible monopolar, Luer-Lock connector, 2,2 mm, 165 cm
	HF22-2301	Suction-koagulation elektrode , flexible monopolar, Luer-Lock connector, 2,2 mm, 230 cm
	SL23-1610	Polypectomy snare , complete with tube, OD:2,2 mm, WL:165 cm, asymmetric snare, 2,5 cm wide, without handle,
	SL23-2310	Polypectomy snare , complete with tube, OD:2,2 mm, WL:220 cm, asymmetric snare, 2,5 cm wide, without handle,
	SL23-1612	Polypectomy snare , complete with tube, OD:2,5 mm, WL:165 cm, symmetric snare, 3,5 cm wide, without handle,
	SL23-2312	Polypectomy snare , complete with tube, OD:2,5 mm, WL:220 cm, symmetric snare, 3,5 cm wide, without handle,
	HS11-1200	Handle for Polypectomy snare , standard size
	HF-301-ER	Standard HF-cable , monopolar, 3 m

Biopsy/grasping forceps

	article no.		Ø (mm)	length (cm)
	B22-16011	Biopsy forceps with spike	2,2	160
	B22-22011	Biopsy forceps with spike	2,2	220
	F22-16041	Grasping forceps	2,2	160
	F22-22041	Grasping forceps	2,2	220
	F22-16040	Grasping forceps	2,2	160
	F22-22040	Grasping forceps	2,2	220
	F23-16042	Grasping forceps	2,3	160
	F23-22042	Grasping forceps	2,3	220
	FK18-1635	Basket, 4 wires	1,8/3,5x5cm	165
	FK23-1635	Basket, 6 wires	2,2/4,5x6cm	165
	FK23-1637	Basket, 6 wires	2,2/5x10cm	165
	FK23-2035	Basket, 6 wires	2,3/5x6cm	220
	X-46-00	Foreign body protection cap for gastrofiberscopes, Ø 9-11 mm		