

Dr. Fritz®

GmbH

Endoscopes and Video Systems

Vet-Info 12-1 E

Dogs

Equipment for gynaecological and urological examination in small animals acc. to Blendinger



Canine gynaecology – Cystoscopy



Endoscopes for obstetrical work in small animals

Vaginoscopy is a common procedure in small animal practices. However, there are practical difficulties in small bitches or due to pain reaction depending on the animals cycle. The new "Blending" video vaginoscope set is the solution for gynaecological examination of all sizes of bitches. Small dogs will tolerate the 5 mm outer diameter of the examination sheath, larger instrumentation is resented. The DIN-ocular telescope is compatible to all standard endoscopic video cameras. There is also a so called Twin connector telescope available which fits

directly to the new video compact system. The Twin connector cable carries both light and electronic signal reducing clutter and tangles.

In addition to, and a consequence of, being patient friendly, the resulting on screen image allows the surgeon to fully explain to the owner exactly what is happening. If the owner understands what is required compliance (and payment!) is much more likely. With the newly developed working sheath artificial insemination in dogs has reached a new dimension. The scope and sheaths are also useful for cystoscopy and other endoscopic procedures.



- | | |
|-------------------|--|
| T301-4030 | Wide angle telescope, 30°, OD:4mm, WL:30cm, 115° field of vision, autoclavable, standard light cable connectors |
| T301-4030W | Wide angle telescope, 30°, OD: 4mm
Twin Connector, 30cm, 115° field of vision, comp. to video compact system |
| S301-4034 | "Blending" vaginoscope sheath I,
for video-vaginoscopy and urology, OD:5mm, WL:27cm,
comp. with telescope T301-4030/W, inflation stopcock |



- | | |
|------------------|---|
| S301-4037 | "Blending" vaginoscope sheath II,
17 Chrr/Fr working sheath for vaginoscopy and urology, WL:25 cm
comp. with T301-4030/W, inflation stopcock, 2mm working channel. |
|------------------|---|

Accessories:



- | | |
|------------------|--|
| H150-1-1 | Inflation pump, small, 5 litre/min, max. 0,25 bar,
230 VAC, 50/60Hz, incl. silicon tube, Luer Lock |
| F11-69 | Special catheter 4 Fr., 69 cm, olive end, angled tip, y-sterile |
| F12-69 | Special catheter 5 Fr., 69 cm, olive end, angled tip, y-sterile |
| F13-69 | Special catheter 6 Fr., 69 cm, olive end, angled tip, y-sterile |
| F16-50 | Special catheter OD:1,6mm, 50cm
autoclavable, straight, rounded end. |
| B18-33 | Flexible biopsy forceps, OD:1,8mm,
35 cm comp. to rhinoscope and vaginoscope |
| TR-60-1 | Sealing cap. for small trocar cannulas
cap opening D:1mm, grey, 5pcs |
| K-3325-62 | Transportation case |