




Endoscopes and Video Systems

Vet-Info 21-4 E

Suction-/irrigation pump the perfect combination for the small and large animal practice



Suction-/irrigation pump

the perfect combination for the small and large animal practice

- *suction and irrigation separately or at the same time*
- *individually adjustable*



H160-2S Veterinary suction & irrigation pump, complete with 2 bottles (suction/ water) 1 litre each, silicon tubing set, vacuum: 11 l/min, 0-08bar, pressure: 0-0,7bar (up to 140mmHg) fully adjustable, 230VAC;50/60Hz, incl. filters



AR-10-922 Spare filters, sterile, for CO² infusor and air pressure, disposable, box with 50 pcs.