

Dr. Fritz[®]

LLC

Endoscopes and Video Systems

Vet-Info 18-1 E

LED
Dentalstick
Equine

Equine, Bovine, Camels, Large Animals

Equipment for endoscopic examination of teeth
in field practice
– wireless option –



wireless




**LED Dentalstick Equine
for large animals – digital**

LED Dentalstick Equine

the revolution in veterinary dental care

- the scope for equine dentists
- wireless technique ideal for field practice
- easy handling, quick cleaning and disinfection
- perfect presentation with the owners
- video transmission wireless or via USB/video cable
- 1/6" CCD chip
- Resolution 750x580 pixel



Basic System with video and USB connection



Monitor 5,6"

wireless



Sender/Battery Unit

Equine dental camera with LED light – basic system

TW61-150B Dentalstick-Equine "Compact", dental video camera with integrated LED light, field of vision: 115°, OD:15mm, WL:52cm, connections: video/USB out, including video cable, power transformer (115-230VAC), software for PC/notebook
transportation case - compatible to wireless modus

Portable monitor, small, analog

M-56-A LCD-Monitor 5,6", sunlight readable, power transformer, video/cable, 0,8 kg

Equine dental camera with LED light and wireless video transmission (<15m)

TW61-150S Dentalstick-Equine "Compact-W", SET with wireless option

consisting of:
Dentalstick Compact (TW61-150B)
Wireless Modul Set incl. sender and USB receiver unit (TW61-150W)