




Endoscopes and Video Systems

Vet-Info 17-6 E


Equine, Bovine, Camels, Large Animals

Flexible gastroscopes for examination of the respiratory system and the GI tract

Gastroscopy
Bronchoscopy



Bronchoscopy/Gastroscopy




Gastroscopy/Bronchoscopy for large animals



Video-Bronchoscopes and Gastrosopes



Completely mounted



- FVE-130-10** **Equine video bronchoscope, 1/6" CCD chip,** OD:9,8 mm, WL:135 cm, full size screen technology, standard buttons, freeze function, angle of view 115°, working channel 2,8 mm, standard connectors for cleaning and disinfection, light connector, compatible to Olympus, including density tester, cleaning brush, processing cable, case.
- FVE-190-10** **Equine video bronchoscope, 1/6" CCD chip,** OD:11,5 mm, WL:185 cm, full size screen technology, standard buttons, freeze function, angle of view 115°, working channel 2,8 mm, standard connectors for cleaning and disinfection, light connector, compatible to Olympus, fully immersible, including density tester, cleaning brush, processing cable, case.
- FVE-300-13** **Video gastroscope for large animals, 3 m, 14" CCD** same but OD: 12,9 mm, WL: 330 cm

4 way deflection

Video-Processor, Light source

– perfect illumination for gastroscopy, laparoscopy and arthroscopy –

Vet Video Combo Xe 100,

– universal combination of processing unit and xenon light –

- easy to use, for all applications
- economical

48-P28-K



48-PXE-100


- 48-PXE-100** **Vet Video Combo Xe 100,** digital CCU & 100 w Xenon light source, for FVE video gastro/colonoscopes and rigid video cameras, PAL, 110-240VAC, 50/60Hz
 - light connector comp. to Olympus
 - built-in air pump
 - automatic exposure control, 2mill.shutter
 - 2,5 Lux min, light sensitivity
 - resolution 752(H)x582(V), 440 000 pixel
 - video outlet: VHS, S-VHS,

Accessories, optional:

- 48-P28-K** Camera head
- 48-XE-100E** Xenon spare lamp type XBO R

Suction-/irrigation pump

the perfect combination for the small and large animal practice



H160-2S

Veterinary suction & irrigation pump, complete with 2 bottles (suction/ water) 1 litre each, silicon tubing set, vacuum: 11 l/min, 0-0,8bar, pressure: 0-0,7bar (up to 140mmHg) fully adjustable, 230VAC;50/60Hz, incl. filters

AR-10-922

Spare filters, sterile, for CO² infusor and air pressure, disposable, box with 50 pcs.



- suction and irrigation separately or at the same time
- individually adjustable



WV-53022

Video cart, size:62x63cm, h:110cm, completely mounted, basic version, color:grey, partly blue,

- drawer inbuilt
- specially rust protected (without holder for gastroscopie and infusion)

Accessories:



47-DK-602

Connecting cord, 18 pin, Dr-Fritz CCU to Dr.Fritz video gastro/colonoscopes, 40cm

X-10-27

Irrigation bottle, 250ml, for video bronchoscope and gastroscopes

EF-18

Density tester comp.to flexible video gastro / colonoscopes, Typ FVE

RH99-56

Cap for instrument channel, bag with 5 pcs

F28-1709

Cleaning brush, flexible, 160cm for working channels 2,6- 3,5mm

F28-2208

same but with 220cm

PC-95-PSK

USB kit, software driver with cable connection, Chinch/Y/C-USB

Accessories:

Catheter, injection needle



- F20-220** Catheter f. gastro fiberscopes, Teflon, OD:2mm, WL:235cm, LL(f) connector, dismantable, autoclavable
- F20-219** Spare catheter f. gastro fiberscopes, Teflon, OD:2mm, 10 m, autoclavable
- LLW-1220** Luer-Lock-connector f. tube, 2mm, luer-lock female,




- F231-220** Flex. injection needle, OD:2,2 mm, 240 cm, Teflon tube, needle size: 0,6x10mm, LL(f)


Accessories for high frequency surgery




- HF22-2301** Suction-koagulation elektrode, flexible monopolar, Luer-Lock conector, 2,2 mm, 230 cm



- SL23-2310** Polypectomy snare, complete with tube, OD:2,2 mm, WL:230 cm, asymmetric snare, 2,5 cm wide, without handle,




- SL23-2312** Polypectomy snare, complete with tube, OD:2,2 mm, WL:230 cm, asymmetric snare, 1,5 cm wide, without handle,




- HS11-1200** Handle for Polypectomy snare, standard size

- HF-301-ER** Standard HF-cable, monopolar, 3 m

Biopsy/grasping forceps



article no.	Ø (mm)	length (cm)
B22-16011 Biopsy forceps	2,2	160
B22-21011 with spike	2,2	230
B23-3601		2,3 360



F22-16041 Grasping	2,2	160
F22-22041 forceps	2,2	230
F23-3604		2,3 360



F22-16040 Grasping	2,2	160
F22-22040 forceps	2,2	220



F23-16042 Grasping	2,3	160
F23-22042 forceps	2,3	230



FK23-1635 Basket, 6 wires	2,2/4,5x6cm	160
FK23-1637 Basket, 6 wires	2,2/5x10cm	160
FK23-2035 Basket, 6 wires	2,3/4,5x6cm	230