


Dr. Fritz®

Gas-4-01E

Flexible LED gastroscopes/coloscopes

for examination of
the GI tract



- Digital LED technology
- Battery operated
- Wireless video transmission
- Easy to use and uncomplicated handling

only for export outside EU
– non registered product –


Wireless LED Video Gastrosopes/Coloscopes

Switch on scope and receive image on a monitor or PC, wireless
Cleaning and disinfection is much easier due to less connection parts

- video gastroscope with integrated and LED-Light
- resolution 750x580 pixel
- integrated sender for wireless video transmission
- easy handling, quick cleaning and disinfection
- perfect video presentation



– possible connections –



- TXH106-9E Wireless LED video gastroscope**, WL: 105cm, OD: 9,8mm, in-built video camera with LED
working channel: 2,8mm, 4-way deflecting
full screen technology, fully immersible
- in-built 5,4 GHz sender for wireless video transmission to analog monitors and/or notebooks/netbooks/PC's, including:
 - charger | 10-240 VAC
 - density tester
 - cleaning brush
 - cleaning adaptors
 - transportation case
 - biopsy forceps
 - 250ml water bottle and air/water double tube
 - TV receiver (TXLI-581)


TXH141-12E Wireless LED video colonoscope, channel 3,4mm, WL: 135cm, OD: 11mm

TXH161-12E Wireless LED video colonoscope, channel 3,4mm, WL: 160cm, OD: 11mm

only for export outside EU
– non registered product –

Optional accessories:

Video Connections – wireless –



Receiver:



TXLI-581

TXLI-7-NO
Receiver Notebook/PC

PC-95-4EU

TXLI-581 TV receiver, 5,8GHz, only, standard, for TXL/TXE series, video outlet: VHS/FBas, 60g, power supply: 115-240VAC/50/60Hz

TXLI-7-NO Receiver, 5,8GHz for notebooks and PC s, with power supply,

***installation by the customer

PC-95-4EU Recorder notebook, HD Endostore software, completely installed incl. database and video cut, compatible "ready to go", compatibly to network systems, Win 7, with powersupply and EU-Keyboard

Spare parts:



EF-18 Density tester
comp. to flexible video gastro- /colonoscopes,



RH99-56 Cap for instrument channel, bag with 5 pcs



F28-1709 Cleaning brush, flexible, 160cm for working channels 2,2-3mm



TW4-250 Air/water bottle, 250ml for LED video-broncho/gastrosopes, autoclavable




TW4-133 Air/water twin tube, 2,6m, complete with couplers for TW/TX series

Recommended accessories:

Small air pump

for lens cleaning and inflation




TW4-249

Air pump and holder, for air/water supply in LED gastroscopes, consisting of:
pump TW4-240 and holder TW4-250H

Suction-/irrigation pump

for lens cleaning, inflation and suction

- low noise and soft pressure system
- easy handling
- separate combined use of suction and air pressure



HI 40-3SH

Suction & irrigation pump, soft system, complete with 2 bottles (suction/ water) 1 litre each, tubing set autoclavable, vacuum: 3 L / min, max: 0,4 bar pressure: 4 L / min, max.0,5 bar (400mmHg) incl. transformer 115- 230VAC; 50/60Hz

AR-10-922

Sterile filter for CO² infuser and air with ID: 4-6mm, disposable, 50pcs/box



only for export outside EU – non registered product –